User Manual

Heated Breathing Tube

LH1

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


CE 0123

Table of Contents

1. Device Symbols	1
2. Warning, Caution and Important Tip	2
3. Intended Use	2
4. Intended Clinical Benefits	2
5. Specifications	3
6. Features	4
7. First Time Setup	4
7.1 Connecting the LH1 Tube to the Device	5
7.2 Removing the LH1 Tube from the Device	6
7.3 Requirements for breathing gas supply device	6
8. Cleaning and Maintenance	6
8.1 Cleaning the LH1 Tube	6
8.2 Disinfection of the LH1 Tube	7
9. Disposal	7
10. Technical Support	8
11. EMC Requirements	8
12. Limited Warranty	12

1. Device Symbols

	Caution		Logo of BMC Medical Co., Ltd.
	Non-sterile		CE marking
	Do not use if package is damaged and consult instructions for use		Made in China
	Importer		Manufacturer
	Batch code		Use-by date
	Model number		Temperature limit
	Medical device		Humidity limitation
	Unique device identifier		Atmospheric pressure limitation
	Authorized Representative in the European Community		DC Power
	Follow instructions for use		Type BF Applied Part
	Dripping (15° tilted)		Marking of electrical and electronic equipment

2. Warning, Caution and Important Tip

WARNING!

Indicate the possibility of injury to the user or operator.

CAUTION!

Indicate the possibility of damage to the device.

IMPORTANT TIP!

Indicate the possibility that such operation may affect the effectiveness or ease of use of the device.

Warnings, Cautions, and Important Tips appear throughout this manual as they apply.

3. Intended Use

The LH1 Heated Breathing Tube (LH1 Tube for short) is heated air tubing intended for incorporation into respiratory therapy devices and humidifiers with conical ISO connectors (ISO 5356), and is intended to reduce or eliminate water condensation and / or pooling of water in the LH1 Tube, and problems associated with such.

The LH1 Tube is indicated for non-invasive respiratory therapy in the home, hospital or sleep-lab setting by a single adult patient. It can also be used in conjunction with supplemental Oxygen.

WARNINGS!

- This device is intended for adult use only.
- This product is not intended for patients with upper respiratory tract bypass.
- The device is intended for use by a single patient only.
- The instructions in this manual are not intended to supersede established medical protocols.
- When the product is connected with other devices to constitute the ME system, it shall meet the relevant requirements of IEC 60601-1-1 for the safety of medical electrical system.
- This device has no alarm function and is not applicable to the case where the alarm is required.
- Images shown in this manual are indicative only. If there is inconsistency between the image and actual product, the actual product shall govern.
- Any serious incidents that occur in relation to this product should be reported to BMC and the competent authority in your country.

4. Intended Clinical Benefits

- (1) Improve patients' quality of sleep.
- (2) Improve patients' compliance with therapy.

5. Specifications

Inner Diameter	19 mm		
Length	The length under static (unstressed) condition in the horizontal plane is 1.8 m Margin of Error: $\pm 10\%$		
Environmental Conditions		Operation	Transport and Storage
	Temperature	5°C to 35°C	-25°C to 70°C
	Humidity	15% to 93%, Non-condensing	15% to 93%, Non-condensing
	Atmospheric Pressure	760 hPa to 1060 hPa	760 hPa to 1060 hPa
Power Supply	DC Voltage		Load Power
	24 V		18 W
Air outlet Temperature	At temperature 5°C to 30°C, air outlet temperature is less than 41°C.		
Resistance to flow at rated flow	30 L/min: < 0.06 hPa/L/min		
Compliance standards	60 hPa: < 5 mL/hPa		
Manufacture Date	Refer to the nameplate		
Service Life	The product's service life is six months from the first day of use, if the use, maintenance, cleaning and disinfection are in strict accordance with the User Manual.		
Validity Term	3 years		
Type of Protection against Electric Shock	Not Applicable		
Degree of Protection against Electric Shock	Type BF Applied Part		
Degree of Protection against Ingress of Water	IPX2		

Note: The nominal rated flow of the product is 60 L/min, 0.2 kPa.

6. Features

The LH1 Tube is composed of breathing gas supply device connector, mask connector, power connector, threaded tubing.

The structure of the LH1 Tube is shown in Fig. 1.


Fig. 1

WARNING!

- This product should only be connected with non-invasive breathing gas supply device or positive pressure breathing gas supply device (hereinafter referred to as breathing gas supply device) produced or authorized for use by BMC or recommended by prescription doctor. The use of incorrect breathing gas supply device may affect the the effectiveness of the treatment.

7. First Time Setup

WARNINGS!

- When the indoor temperature exceeds 30°C, disconnect power supply of the LH1 Tube, so as to avoid the high temperature of airflow and irritate the airway of patients.
- When the indoor humidity exceeds 80%, do not use the LH1 Tube to prevent any condensation in the tubing from entering the patient's airway.
- Do not cover the LH1 Tube with anything, including textile heat insulation sleeves, plastic sleeves, blankets, etc., to avoid excessive air temperature and irritate the airway of patients.
- Do not use the LH1 Tube without airflow.
- When accidents occur during normal use, the product should be disconnected immediately and appropriate emergency and corrective measures should be taken.
- If the DC voltage exceeds the range (refer to Part 5 "Power Supply"), the product will not work properly.
- If the LH1 Tube is damaged (broken hole, kink, tear, exposed heated wire, etc.) or poor function, please do not repair and use it by yourself and replace it immediately.

CAUTIONS!

- Check the integrity of the product packaging, if the package is found damaged, please do not use this product.
- When using the product, make sure that patients are closely monitored.
- Check the product regularly and replace the product immediately if it is found to be contaminated or no longer applicable to its use intention.
- Do not pile too long tubs on the head of the bed, it is possible to entangle the head or neck of the patient in a deep sleep. The tubing should not be covered by the sheet or affected by the heating source (such as electric blanket), otherwise it may be deformed and cause danger.
- When the product is not connected to the patient, avoid foreign matter entering the inside of the LH1 Tube.
- Before each use, the tubing should be checked for damage or foreign matter. If so, please replace the tubing immediately. Check whether the mask and the tube system is sealed.
- When the product is discarded, the treatment method of the high-molecular medical device shall be carried out for the harmless treatment.
- When connecting or disconnecting LH1 Tube with breathing gas supply device or nasal cannula, be sure to hold on the adapter and do not pull the threaded tubing.
- When using the LH1 Tube, be sure to use with the humidifier, so as to avoid the gas temperature being too high.
- When using electrical products, be sure to follow the usual safety measures.
- Do not use (microwave) oven, fan or any other appliances to dry the cleaned LH1 Tube. Please refer to the cleaning instructions.
- Do not use a brush or any other object to clean the inner wall of the LH1 Tube.
- Do not disassemble or intentionally damage the LH1 Tube.
- When you want to permanently discard the LH1 Tube, please clearly mark the used LH1 Tube to avoid misuse by other patients.
- Do not use this product in an environment with flammable gases (such as narcotics) to avoid explosion.

7.1 Connecting the LH1 Tube to the Device

- (1) Unpack and take out the product.
- (2) Connect the breathing gas supply device connector of LH1 Tube and breathing gas supply device.
- (3) Connect the power connector of LH1 Tube and the socket of breathing gas supply device as shown in Fig. 2.


Fig. 2

- (4) Connect the mask connector of LH1 Tube and mask.

CAUTION!

- Before each use of the product, test the product according to the instructions of the device manufacturer to ensure that the product is compatible with the device and that the product is free of air leakage and blockage.

7.2 Removing the LH1 Tube from the Device

- (1) Close the breathing gas supply device.
- (2) Press and hold the power connector of the LH1 Tube to separate it from the socket on the breathing gas supply device.
- (3) Disconnect the LH1 Tube from breathing gas supply device and accessories.

7.3 Requirements for breathing gas supply device

- (1) The breathing gas supply device should have current limitation when supplying power to the LH1 Tube, and the output power supply should meet the requirements of "Power Supply" in part 5 of this manual.
- (2) The socket on the breathing gas supply device must match the power connector of the LH1 Tube.

8. Cleaning and Maintenance**WARNINGS!**

- Regular cleaning of the LH1 Tube is very important for the prevention of respiratory infections.
- To avoid electric shock, always unplug the device before cleaning.
- Use mild soap that is nontoxic to humans.
- Do not use irritating soap, detergent, solvent or alcohol clean the tubing.
- Do not reconnect the LH1 Tube until it is completely dry.
- Do not wash LH1 Tube with dish-washing machine.
- The power connector of the LH1 Tube should not have any contact with any liquid and a disposable medical cleaning cloth or paper is used to remove any visible dirt.
- The LH1 Tube can be cleaned up to 180 times and disinfected up to 100 times.

8.1 Cleaning the LH1 Tube

Preparation before cleaning

- (1) Hold the tube connector and disconnect the tube from the device.
- (2) Prepare tools: soft bristle brush, sterile water (five litres per tube), mild liquid detergent (for example: dilute the soap with clean water to make 2% soapy water).

Instructions for cleaning

- (1) Soak the LH1 Tube in 2% soapy water at 35°C to 40°C for 10 minutes.
- (2) Then clean the inside and outside surfaces with a soft bristle brush 5 times.
- (3) Finally, rinse the inside and outside surfaces with sterile water 3 times, 1 minute each time.
- (4) Gently tap the tube to remove excess moisture from the connector ports. Hang the LH1 Tube so that both ends of the opening face the ground. Let the LH1 Tube naturally dry, out of direct sunlight.

WARNINGS!

- Please wash by hand.
- If the LH1 Tube is damaged (broken hole, tear, exposed heated wire, etc.) or poor function, please do not repair and use it by yourself and replace it immediately.
- Failure to clean in accordance with the Manual may result in reduced performance of the LH1 Tube or reduced product life.

Note: The tube is recommended to be cleaned once a day or after each use.

8.2 Disinfection of the LH1 Tube

If the cleaning instructions in Clause 8.1 are strictly followed, it is generally not necessary to disinfect the LH1 Tube. When the product is contaminated, please clean it according to the cleaning method in Clause 8.1, before disinfecting according to the following scheme.

- (1) Immerse the LH1 Tube in 0.55% O-phthalaldehyde for 12 min.
- (2) After immersing, rinse it in clean water thoroughly.
- (3) Tap gently to remove excess moisture from the connector port.
- (4) Suspend the LH1 Tube with the openings at both ends facing the floor.
- (5) Air-dry the tubing, and avoid direct sunlight.

WARNING!

- No sterilization is required by the product.

CAUTIONS!

- Disinfectants generally damage the surface of materials and shorten the service life of the product. Attention should be paid to the instructions for use of the disinfectants; the right ones should be selected as applicable as possible, and the recommendations of the disinfectant manufacturer should be followed.
- If the LH1 Tube is found to have any damage (such as fracture and rupture) after disinfection, it needs to be replaced in time. It is normal for the LH1 Tube to become slightly discolored.
- After disinfection, please be sure to always rinse it thoroughly with clear water to prevent residual disinfectant from damaging skin and respiratory tractor causing allergic reaction.

9. Disposal

This product complies with the requirements of Directives 2011/65/EU and (EU)2015/863(RoHs) on the restriction of hazardous substances.

According to the European Directive 2012/19/EU (WEEE) on electrical and electronic equipment, products cannot be discarded indiscriminately and must be disposed of in accordance with the laws and regulations of the country in which disposal occurs.

The crossed-out wheeled bin symbol  indicates that the product bearing this symbol may not be disposed of together with general household waste, but instead requires separate disposal. This requirement for separate disposal is based on the European Directive 2012/19/EU for electrical and electronic equipment. You can hand in the product at a municipal collection point, for example. This reduces the impact on natural resources and prevents contamination of the environment through the release of hazardous substances.

For further information regarding product disposal, please contact your specialist distributor.

10. Technical Support

Please contact BMC directly if you need the technical documents of the LH1 Tube for certain purposes such as maintenance or connection to other equipment. BMC will provide the technical documents in whole or in part according to your needs.

11. EMC Requirements

Guidance and manufacturer's declaration - electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The user of the device should ensure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes
Harmonic emissions IEC 61000-3-2	N/A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	N/A	

Recommended separation distances between portable and mobile RF communications equipment and the device			
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output of transmitter (W)	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.80	3.80	7.28
100	12	12	23
Note 1: At 80 MHz and 800 MHz, the higher frequency range applied. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			

Guidance and manufacturer's declaration - electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The user of the device should make sure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floor should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient / burst IEC 61000-4-4	±2 kV for power supply lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment
Surge IEC 61000-4-5	±1 kV line(s) to line(s)	±1 kV line(s) to line(s)	Mains power quality should be that of a typical commercial or hospital environment
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% U_T ; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U_T ; 1 cycle 70% U_T ; 25 / 30 cycles At 0° 0% U_T ; 250 / 300 cycles	0% U_T ; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U_T ; 1 cycle 70% U_T ; 25 / 30 cycles At 0° 0% U_T ; 250 / 300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery
Power frequency (50 Hz / 60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment
Note: U_T is the AC mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration - electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The user of the device should make sure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 V (Effective value) 150 kHz to 80 MHz	3 V (Effective value)	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitter, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	
Note 1: At 80 MHz and 800 MHz, the higher frequency range applied. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device. ^b Over the frequency range 150 kHz to 80 MHz, the field strengths should be less than 3 V/m.			

WARNINGS!

- During operation of the device, due to electrostatic interference, the following phenomena may occur:
 - (1) Temporary loss of function or degradation of performance, such as abnormal screen display, etc. The device will recover to normal after being restarted.
 - (2) Automatic restart of the device. These phenomena will not affect the normal use of the device, and will not cause permanent performance degradation or function loss of the device.
- This device should not be used in the vicinity or on the top of other electronic equipment such as cell phone, transceiver or radio control products. If you have to do so, the device should be observed to verify normal operation.
- The use of accessories and power cord other than those specified, with the exception of cables sold by the manufacturer of the equipment or system as replacement parts for internal components, may result in increased emissions or decreased immunity of the equipment or system.
- Even if other equipment meets the emission requirements of corresponding national standards, the product may still be interfered by other equipment.
- This product is not intended for use with MRI equipment and is not allowed to be used in MR environments.
- Do not use the device or accessories in an environment with electromagnetic equipment such as CT scanners, Diathermy, RFID and electromagnetic security systems (metal detectors) as it may cause unacceptable risk to the patient or damage to the device. Some electromagnetic sources may not be apparent, if you notice any unexplained changes in the performance of this device, if it is making unusual or harsh sounds, disconnect the power cord and discontinue use. Contact your home care provider.
- When the product is exposed to soldering, electrosurgery, defibrillation, X-ray (γ ray), infrared radiation, transient electromagnetic field, including nuclear magnetic resonance (MRI) and radio interference environment, the product may be damaged.

12. Limited Warranty

BMC Medical Co., Ltd. warrants that the device shall be free from defects of workmanship and materials and will perform in accordance with the product specifications for three (3) months from the date of sale by BMC Medical Co., Ltd. to the dealer. If the product fails to perform in accordance with the product specifications, BMC Medical Co., Ltd. will repair or replace, at its option, the defective material or part. BMC Medical Co., Ltd. will pay customary freight charges from BMC Medical Co., Ltd. to the dealer location only. This warranty does not cover damage caused by accident, misuse, abuse, alteration and other defects not related to material or workmanship.

BMC MEDICAL CO., LTD. DISCLAIMS ALL LIABILITY FOR ECONOMIC LOSS, LOSS OF PROFITS, OVERHEAD OR CONSEQUENTIAL DAMAGES WHICH MAY BE CLAIMED TO ARISE FROM ANY SALE OR USE OF THIS PRODUCT. SOME STATES DO NOT ALLOW THE EXCLUSION OR LIMITATION OF INCIDENTAL OR CONSEQUENTIAL DAMAGES, SO THE ABOVE LIMITATION OR EXCLUSION MAY NOT APPLY TO YOU.

To exercise the rights under this warranty, contact the local authorized dealers or:

MANUFACTURER:

BMC Medical Co., Ltd.

Room 110 Tower A Fengyu Building, No. 115 Fucheng Road, Haidian, 100036

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Shanghai International Holding Corp. GmbH (Europe)